




3DHISTECH

# **RENEWED DIGITAL TISSUE MICROARRAY (TMA) MACHINERY AND SOFTWARE WITH PCR EXTRACTION FUNCTION**



**Attila Szautner**




# TMA technology

3DHISTECH

Glass slide

Marked area

1.






# TMA technology

3DHISTECH

Glass slide


Marked area

1.



Matching  
area on the  
donor block

2.






# TMA technology

3DHISTECH

Glass slide


Marked area

1.




Matching  
area on the  
donor block

2.



Embedding

3.




Recipient block



# TMA technology


3DHISTECH





# TMA technology

3DHISTECH





3DHISTECH

# Advantages

- **Save:**

- Reagent
- Time
- Tissue
- Workforce



- **Standardized laboratory conditions**
- **Central data handling**





# Case Selection

3DHISTECH

- Collecting donor blocks
- Preparing clinical data

	A	B	C	D	E
1	<b>Donor Block ID</b>	<b>Organ</b>	<b>Diagnosis</b>	<b>Grade</b>	<b>Sender institute</b>
2	1023-03	breast	IDC	G1	Bay
3	1199-02	breast	IDC	G3	S.E. I. Pathology
4	1222-03	breast	IDC	G2	Bay
5	1315-02	breast	IDC	G2	S.E. I. Pathology
6	1324-03	breast	ILC	G1	Bay
7	1416-04	breast	ILC	G2	Bay
8	1610-04	breast	IDC, comedo	G2	Bay
9	1643-04	breast	IDC	G3	Bay
10	1722-02	liver	cirrhosis	G0	Szombathely
11	1826-02	breast	DCIS	G0	S.E. I. Pathology
12	1877-02	breast	ILC	G2	Natl. Inst. Cancer
13	210-04	breast	Mucinosus carcir	G2	Bay
14	28-05	lymph node	ADC metastasis	poor	Bay
15	322-02	breast	IDC	G2	S.E. I. Pathology
16	33-05	breast	IDC	G3	Bay
17	4761-02	breast	DCIS	G0	S.E. I. Pathology








3DHISTECH

# Analysis of digital slides, Creating markers





3DHISTECH

# Workflow, Evaluation of TMA slides

## EVALUATION & ANALYSIS

### Case Center



### SCANNING THE TMA SLIDES



### Slide preparation tools




### DNA isolation, DNA amplification, DNA sequencing



### PCR extraction

### TMA BLOCK CREATION





# TMA Module Software

3DHISTECH

- **The TMA Module deals with simultaneous handling of many slides and cores**
- **Multi-user option (Consolidator function)**
- **Filter the spots in the gallery by different parameters (unscored – rescored)**
- **Statistical analysis**
- **Simple or Combined scoring scheme**






3DHISTECH

# Workflow steps of the TMA Module

- **Preparing the Excel sheet for the TMA project**
- **Spot finding**
- **Scoring and managing slides**
- **Export data and images**





# Preparing the Excel sheet for TMA project

3DHISTECH

	A	B	C	D
1	Unique TMA Sample ID	TMA Block ID	Donor Block ID	
2		1 Block1	xyz1	
3		2 Block1	xyz2	
4		3 Block1	xyz3	
5		4 Block1	xyz4	
6		5 Block1	xyz5	
7		6 Block1	xyz6	
8		7 Block1	xyz7	
9		8 Block1	xyz8	
10		9 Block1	xyz9	
11		10 Block1	xyz10	
12		11 Block1	xyz11	
13		12 Block1	xyz12	
14		13 Block1	xyz13	
15		14 Block1	xyz14	
16		15 Block1	xyz15	
17		16 Block1	xyz16	
18		17 Block1	xyz17	
19		18 Block1	xyz18	
20		19 Block1	xyz19	
21		20 Block1	xyz20	
22		21 Block1	xyz21	
23		22 Block1	xyz22	
24		23 Block1	xyz23	
25		24 Block1	xyz24	
26		25 Block2	xyz25	
27		26 Block2	xyz26	

	A	B	C	D	E	F
1	Block1	X1	X2	X3	X4	
2	Y1	1	2	3	4	
3	Y2	5	6	7	8	
4	Y3	9	10	11	12	
5	Y4	13	14	15	16	
6	Y5	17	18	19	20	
7	Y6	21	22	23	24	

Block1

	A	B	C	D	E	F
1	Block2	X1	X2	X3	X4	
2	Y1	25	26	27	28	
3	Y2	29	30	31	32	
4	Y3	33	34	35	36	
5	Y4	37	38	39	40	
6	Y5	41	42	43	44	
7	Y6	45	46	47	48	

Block2

	A	B	C	D	E	F
1	Block3	X1	X2	X3	X4	
2	Y1	49	50	51	52	
3	Y2	53	54	55	56	
4	Y3	57	58	59	60	
5	Y4	61	62	63	64	
6	Y5	65	66	67	68	
7	Y6	69	70	71	72	

Block3

	A	B	C	D
1	Stain	Block1	Block2	Block3
2	Ki67	23-Ki67-Mib1	24-PKi67	
3	ER	21-ER-6F11	27-ER-6F11	30-ER-SP1-Chem



# Scoring and managing slides

- Combined Scoring

Scoring Scheme Setup Window

Stains

Selected Stain: K167

Scoring Scheme Type

Simple

Combined

Save Scheme

Number of Columns: 4

Number of Rows: 3

Title of Columns: Frequency

Title of Rows: Intensity


	-10%	10-40%	40-80%	80%--	
Weak	1	2	3	4	+1
Moderate	5	6	7	8	+2
Strong	9	10	11	12	+3
	+1	+2	+3	+4	



3DHISTECH

# Scoring and managing slides, Gallery

- Scoring spots
- Ordering spots
- Filtering spots





3DHISTECH

# TMA Grand Master

- **72 block capacity:**
  - 12 recipient blocks, 60 donor blocks
- **4 donor punch and drill bit sizes:**
  - 0.6, 1, 1.5, 2 mm**
- **2 digital cameras (1 for label and 1 for block view)**
- **Automated block and tool height measurement by sensor**
- **Separately controlled puncher and drill bit**
- **Parallel block feeding, drilling, punching and image capture**







3DHISTECH

# TMA Grand Master

- **Faster movements and faster workflow**
- **Computer-aided tasks:**
  - Recipient block layout design
  - Assign a digitized slide to a donor block
  - Data export to Excel





3DHISTECH

# TMA Grand Master






3DHISTECH

## Jpg image overlay


- **Set the transparency near the block image**





# Barcode reading

- **Auto Barcode button**
- **1D Barcode Types**
  - Code 128, Code 93, EAN8, EAN13, UPC\_A
- **2D Barcode Types**
  - Datamatrix, Maxicode, Micropdf417, QRcode






# PCR in Grand Master

3DHISTECH

- **Dedicated position**





# Cleaning block

3DHISTECH

- **Set the cleaning block**
- **Clean before each PCR extraction**

